1. IDENTIFICATION DU FABRICANT

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2. DEVICE IDENTIFICATION AND PRESENTATION OF THE DEVICE IN THE PACKAGE

⁵⁷Co pen point markers are used in association with a gamma-camera in order to mark cue points on an image generated in an examination by gamma camera.

The pen point marker consists of a passivated aluminum rod into which is introduced a source of cobalt 57 sealed in the conical portion (A).

Figure 2 : ⁵⁷Co pen point marker (A = radioactive part)

3. PRESENTATION OF THE RANGE OF PRODUCTS

⁵⁷Co pen point marker is available in four kinds with an activity of 3.7 MBq or 7.4 MBq.

Medical Device reference	Medical Device name
CO57EMSA40	⁵⁷ Co pen point marker, activity 3.7 MBq
CO57EMSA80	⁵⁷ Co pen point marker, activity 7.4 MBq
CO57EMSA40N	⁵⁷ Co pen point marker, nominal activity 3.7 MBq
CO57EMSA80N	⁵⁷ Co pen point marker, nominal activity 7.4 MBq

The Medical Device does not have any associated accessory for the use.

4. INSTRUCTIONS FOR USE

4.1. Intended uses

⁵⁷Co pen point marker is solely used to point to a region of interest on an image obtained by a gamma camera. Its use is restricted to authorized personnel to perform radiology examinations by gamma radiation. The staff shall be empowered to handle sealed sources.

LEA cannot be held liable for an anomaly occurring during use other than that intended.

4.2. Target group and expected patient population

There is not really a target group because the radiologist may decide to use or not (depending on the examination method) the ⁵⁷Co pen point marker in order to point to a region of interest.

4.3. Expected performance

Expected performance of the device is:

- a radioactive activity detected by gamma-camera over a 2 years period * (free of risk for an individual patient),
- the source remains sealed for 10 years (date when all sealed sources must be returned to the manufacturer).

(*) Cobalt 57 activity decreases rapidly over time. The radioactive period is 272 days. After 2 years, the residual activity is only 15% of the initial activity. Thus, in order to get the same result with the gamma-camera, measurement will increase from 6 to 7 times. This is the reason why it is recommended not to use the device after 2 years of use.

4.4. Contraindications

There is no known contra-indication under a gamma-camera examination.

4.5. Side effects and potential complications

There is no known side effect due to the use of a pen point marker during a gamma-camera examination.

4.6. Use precautions

Use precautions are mainly related to the safety of staff and patients.

- ^{57}Co pen point marker must be used with gloves and shall not be placed in direct contact with the patient (the patient is dressed during the examination or the pen point marker is placed on a support),
- Outside the time of use, the pen must be placed in the sealed support constituting the radiological protection and then stored in an appropriate storage,
- Prohibit the patient to manipulate the ^{57}Co pen point marker,
- The shielded portion of the package protects against the gamma ray emission from the source, However, no system is completely safe against gamma rays that's why it is forbidden to keep the marker in the pocket of the lab coat even if the pen point marker is in its sealed packaging,
- If the shielded package is deteriorated, store the ^{57}Co pen point marker in a protective area against gamma rays and please contact your supplier in order to get a new shielded package,
- Keep the calibration certificate at your gamma camera in order to be able to verify that the ^{57}Co pen point marker is still active enough before the medical examination,
- Never attempt to open the sealed source,
- use a conventional disinfectant in order to clean the ^{57}Co pen point marker after each use,
- If you notice a defect performance of the ^{57}Co pen point marker, please contact the manufacturer.

4.7. Lifetime of the device

The certificate which is delivered with the ^{57}Co pen point marker allows to know the activity of the radionuclide on the date indicated on the label (12h UTC, DD/MM/YYYY). Like all radioactive sources, the activity decreases over time according to the half-life period. For ^{57}Co isotope, the radioactive period is 272 days. The activity of the ^{57}Co pen point marker has been determined in order to produce a visible signal on an examination by gamma camera for 2 years after the date indicated on the label. These two years are the recommended duration of use.

Be careful: the sealed source remains radioactive after these 2 years period but does not allow anymore to perform a medical examination under reasonable conditions either for the patient or for the staff.

5. USED MATERIALS

The ^{57}Co pen point marker is made of solid passivated aluminium provided with a screwed and glued end-piece, which contains cobalt 57 salts. The sealing is effected by an epoxy glue.

6. STORAGE AND HANDLING CONDITIONS

Outside the time of use, the ^{57}Co pen point marker must be placed in an appropriate storage like any radioactive source. A biological shielding made of lead is provided in the package for this purpose.

To ensure the cleanliness of your pen, it is important to clean it after each use with a conventional disinfectant.

7. ELIMINATION

After 2 years of use you can either return the sealed source in its original packaging to LEA (Radioactive Standard Laboratory) or keep it with the appropriate securities for a maximum of 10 years after the manufacturing date on the label.

The regulation for sealed sources request that the sources shall be returned to the manufacturer after 10 years. Please note that the packaging must **IMPERATIVELY** be returned with this pen, during the recovery operation after use. The recovery procedure will be sent on request.

8. SYMBOLS USED ON THE DEVICE AND THE PACKAGING (ACCORDING TO NF EN ISO 15223-1)

Symbol used	Description
 <p>LEA Site Orano du Tricastin BP 75 26701 Pierrelatte Cedex - France</p>	Name and address of the manufacturer
	CE marking for European Conformity
	Consult the handling instructions
	Danger radioactive material
 <p>CO57EMSA40 or CO57EMSA80 or CO57EMSA40N or CO57EMSA80N</p>	Catalog reference
 <p>xxxx</p>	Serial number xxxx: unique number of the source
 <p>DATE (12h U.T.C) YYYY/MM/DD</p>	Reference date for the source activity
 <p>YYYY/MM/JJ</p>	Expiry date = date when the source activity < 0,5 MBq

Reference of the handling instructions, revision date and year in which the CE marking was affixed

Reference of the handling instructions: SPE.DMM.21.083/00

Revision date: 22th March 2022 edition. English Version

UE declaration of conformity: SPE.DMM.21.081 dated: 12th August 2021

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